System blue

System components

for Photodynamic Diagnosis (PDD)

System blue

System-Komponenten

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Procedure for Photodynamic Diagnosis (PDD)

This procedure is based on the interaction of light of an appropriate wavelength range with a tumor-selective enriched substance. This substance enables a color and brightness contrast based on fluorescence to be generated in Photodynamic Diagnosis (PDD). Fluorescence involves the property of bodies or materials to convert the light absorbed by them into light with other wavelengths.

A photosensitive tumor marker is required in order to perform PDD for early identification of a urinary bladder carcinoma. The currently authorized marker is dissolved in saline solution prior to cystoscopy, and instilled transurethral into the bladder through a catheter.

The bladder requires approximately 60 minutes to completely absorb the administered substance and convert it into the body's own dye protoporphyrin IX (PPIX).

Tumor cells have a metabolic defect which leads to a selective enrichment of PPIX there. After illumination and excitation with blue-violet light, red or magenta (pink) colored fluorescence light is generated in PPIX.

Tissue that has undergone pathological change or a tumor which is easily overlooked in conventional bladder cystoscopy due to flat growth (carcinoma in situ) or small size – is then visualized in a red fluorescence or magenta colored fluorescence. Healthy tissue mostly appears colorless with a slightly cyan hue.

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- ENDOLIGHT LED blue:** The first PDD LED-based light source offers maximum color contrast for clear tumor differentiation
- Autoclavable Special Fiber Light Cable:** Unique for fluorescence applications
- ENDOCAM Logic HD:** High-resolution camera controller with two different PDD modes – **blue PDD** and **blue PDD Color Contrast**
- PENDUAL blue HD:** Articulated full HD camera head for application in cystoscopy and in resection
- PANOVIEW Telescopes blue:** Improved light transmission and longer life time thanks to new fiber design

ENDOLIGHT LED blue bundle

consisting of:

ENDOLIGHT blue (5165001),
Patch cable RJ45 SFTP, 0.5 m (72325378),
Foot control rocker switch (2030.105)
and power cable (2440.03), 100-240 V~, 50/60 Hz 51650011

also:

Special Fiber Light Cable Set Ø 3.5 mm, 2.3 m long consisting of:
Fiber light cable (80673523), adapter projector side (8095.07) and
adapter endoscope side (809509) 806735231

PANOVIEW telescope blue, Ø 4.0 mm

0°	8650.514
12°	8654.531
30°	8654.522
70°	8650.515

Pendual blue HD Camera Head for ENDOCAM Logic HD
Cable length 3 m 5525833

Note:

Application in **PDD** mode only requires ENDOCAM Logic HD 5525 and software update R19.

The use of a 27" monitor from R. Wolf is recommended.

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